


New! New! New!

The "Gryphon" Liquid Handling System

The Gryphon is a personal sized version of the "Phoenix" liquid handling system. By using a combination of interchangeable multiple syringe heads, pipette tip heads and non-contact dispensers, the competitively priced Gryphon provides increased efficiency and cost-effectiveness. The Gryphon minimizes the amount of reagent volume needed and performs many microtiter plate assays. It is designed to be used in a bench-top lab environment.



Features:

- Interchangeable heads:
96 or 384 flexible needle syringes or **pipette tips**
- Base unit includes touchscreen interface for easy control and integrated wash system in a small 9 x 16 x 25 inch (23 x 41 x 63,5 cm) (W x D x H) foot print.
- Optional X-Y stage includes easy to use Windows based software in a small 21 x 21 inch foot print
- Optional non-contact nano dispenser makes an ideal configuration for **protein crystallization** assays.



Please note also available is the "Phoenix" – the bigger version for high-throughput liquid-handling. Details upon request.